



# QUOTATION

By: Alex Xu

**Item: Shiatsu Rolling Back Cushion**

**Date: 10<sup>th</sup> Apr, 2008**

**Model: DK-163**

**No.: DK080410163**

Function	Picture
<ul style="list-style-type: none"> <li>I Continual speed control</li> <li>I 8 obverse/reverse kneading massage balls</li> <li>I Used for shoulders, back, waist, thighs, calves and feet</li> <li>I Make of wear proof and ventilate cloth</li> <li>I Easy to carry, easy usage at different places</li> </ul>	

Price			
<b>FOB XIAMEN</b>	<b>USD 28.00</b> (Without RoHS)	<b>PAYMENT</b>	By T/T; 30% advance, 70% before shipment.
	<b>USD 31.00</b> (With RoHS)		
<b>Lead Time</b>	35 Days for first order 25 Days for repeat order	<b>MOQ</b>	500 PCS (For No RoHS) 2,400 PCS (For RoHS)
<b>Sample Cost</b>	Twice as the price	<b>Sample Availability</b>	10 Days

**Remark:**

- I **The price includes LOGO printed on product/manual/color box.**
- I **Air freight for samples to be paid by customer.**
- I **USD exchange rate is 7.6-7.7**

Packing Details					
<b>PACKAGE</b>	<b>Color Box Size</b>	44X12.5X46.5 CM	<b>N.W./PC (KGS)</b>	20.0	
	<b>Carton Size</b>	64X47X48 CM/5pcs	<b>GW/CARTON (KGS)</b>	22.0	
	<b>Manual</b>	English Available	<b>FCL</b>	<b>20' Q</b>	1,050 PCS
				<b>40' Q</b>	2,100 PCS
		<b>H. Q</b>		2,400 PCS	
<b>DATA</b>	~220V OR ~110V-120V      Rated power: 48W				
<b>Approval</b>	CE / RoHS (need 60 days to apply RoHS)				